

Instructions for use

## Infinity<sup>®</sup> Acute Care System



### **WARNING**

To properly use accessories, read and comply with all relevant instructions for use.

**Monitoring accessories**  
**Software VG7.n**



## Typographical conventions

---

- 1** Consecutive numbers indicate steps of action, with the numbering restarting with "1" for each new sequence of actions.
- Bullet points indicate individual actions or different options for action.
- Dashes indicate the listing of data, options, or objects.
- (A) Letters in parentheses refer to elements in the related illustration.
- A** Letters in illustrations denote elements referred to in the text.

## Information about this document

---

### Trademarks

---

- Infinity®
  - MPod®
  - MCable®
  - MonoLead®
- are trademarks of Dräger.
- rainbow®
  - Masimo®
- Masimo® SET® (Signal Extraction Technology)
- are registered trademarks of Masimo Corporation.
- Capnoline®
  - Microstream®
- are registered trademarks of Medtronic
- TOFscan®
- is a registered trademark of IDMED

### Definitions

---

#### **WARNING**

**A WARNING statement provides important information about a potentially hazardous situation which, if not avoided, could result in death or serious injury.**

#### **CAUTION**

A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or in damage to the equipment or other property.

#### **NOTE**

A NOTE provides additional information intended to avoid inconveniences during operation.

# Contents

---

<b>Information about this document</b> . . . . .	4
<b>For your safety and that of your patients</b> . . . . .	7
<b>List of accessories</b> . . . . .	11
Infinity Acute Care System accessories overview . . . . .	13
General IACS accessories . . . . .	17
IACS / M540 options. . . . .	21
ECG monitoring – reusable leadsets . . . . .	22
ECG monitoring – single-patient use leadsets . . . . .	24
ECG monitoring – MonoLead . . . . .	26
ECG monitoring – neonatal . . . . .	28
Invasive pressure (IP) with Hemo2 Pod or Hemo4 pod . . . . .	29
IP Monitoring with MPod – QuadHemo or MCable – Dual Hemo. . . . .	30
Non-Invasive blood pressure (NIBP) monitoring – Dräger connectors. . . . .	32
Non-invasive blood pressure (NIBP) monitoring – neonatal . . . . .	34
SpO2 monitoring – Masimo rainbow SET (rainbow sensors) . . . . .	35
SpO2 monitoring – Masimo rainbow SET (RD rainbow SET) . . . . .	37
SpO2 monitoring – Masimo rainbow SET (LNCS) . . . . .	39
SpO2 monitoring – Masimo rainbow SET (RD-SET) . . . . .	42
SpO2 monitoring – Masimo SET (LNCS) . . . . .	45
SpO2 monitoring – Masimo SET (RD SET) . . . . .	48
Masimo SpO2 accessories . . . . .	51
SpO2 monitoring – Nellcor OxiMax . . . . .	52
Cardiac (C.O.) output monitoring . . . . .	54
CO2 monitoring – Mainstream . . . . .	55
CO2 monitoring – Microstream . . . . .	56
Temperature monitoring . . . . .	58
MCable – Analog/Sync . . . . .	62
Scio monitoring . . . . .	63
NMT monitoring . . . . .	65

This page has been left blank intentionally.

## **For your safety and that of your patients**

---

Mandatory reporting of adverse events . . . . .	8
Strictly follow these instructions for use . . . . .	8
General precaution. . . . .	8
Training. . . . .	8
Sterile accessories . . . . .	8
Availability of accessories. . . . .	8
Connected devices. . . . .	9
Installing accessories . . . . .	9

## Mandatory reporting of adverse events

Serious adverse events with this product must be reported to Dräger and the responsible authorities.

## Strictly follow these instructions for use

### WARNING

To maintain patient safety, adhere to all **WARNINGS** and **CAUTIONS** on equipment labels or appearing in the Instructions for use for the Infinity Acute Care System (IACS) applications.

### CAUTION

Read all cleaning instructions (for example, originating from the disinfectant manufacturer and the hospital) carefully *before* cleaning the device. Refer to the chapter entitled “Reprocessing” for device-specific cleaning instructions in the Infinity Acute Care System M540 Instructions for use.

## General precaution

### CAUTION

Even accessories designed to be reused have a limited maximum period of use. Handling and reprocessing can increase wear and markedly shorten maximum period of use (e.g., disinfectant residues can attack the material more intensely during autoclaving). If signs of wear become visible, such as cracks, deformation, discoloration, peeling, etc., affected accessories must be replaced.

## Training

Training for users is available via the Dräger organization responsible (see [www.draeger.com](http://www.draeger.com)).

## Sterile accessories

### CAUTION

Risk of medical device failure and of patient injury.

Do not use sterile-packaged accessories if the packaging has been opened, is damaged or if there are other signs of non-sterility. Single-use accessories must not be reused, reprocessed, or re-sterilized.

## Availability of accessories

Some articles are not available worldwide because they are not approved in all countries.

Technical documentation is available on request.

### WARNING

#### Risk due to incompatible accessories

**Dräger has tested only the compatibility of accessories listed in the current list of accessories. If other, incompatible accessories are used, there is a risk of patient injury due to medical device failure.**

**Dräger recommends that the medical device is only used together with accessories listed in the current list of accessories.**



## Connected devices

### **WARNING**

**Risk of electric shock and of device malfunction**

**Any connected devices or device combinations not complying with the requirements mentioned in these instructions for use can compromise the functional integrity of the medical device and lead to an electric shock. Before operating the medical device, strictly comply with the instructions for use of all connected devices or device combinations.**

## Installing accessories

### **CAUTION**

Risk of device failure

Install accessories to the basic device in accordance with the instructions for use of the basic device. Make sure that there is a safe connection to the basic device.

### **CAUTION**

Risk of incorrect use

The accessories are not available individually. Only one copy of the instructions for use is included in the bulk package and must therefore be kept in a location accessible for users.

This page has been left blank intentionally.


## List of accessories











<b>Infinity Acute Care System accessories overview</b> . . . . .	13	<b>Non-Invasive blood pressure (NIBP) monitoring – Dräger connectors</b> . . . . .	32
<b>General IACS accessories</b> . . . . .	17	NIBP hose . . . . .	32
System components . . . . .	17	Reusable adult/pediatric blood pressure cuffs . . . . .	32
Cables and mouse . . . . .	17	Single-patient use adult/pediatric blood pressure cuffs . . . . .	33
Recorders and printers . . . . .	17	<b>Non-invasive blood pressure (NIBP) monitoring – neonatal</b> . . . . .	34
Power supplies . . . . .	17	NIBP hose . . . . .	34
Keyboards . . . . .	18	Single-patient use neonatal blood pressure cuffs . . . . .	34
Network panels and switches . . . . .	18	<b>SpO2 monitoring – Masimo rainbow SET (rainbow sensors)</b> . . . . .	35
Device interface kits . . . . .	19	Masimo rainbow SET intermediate cables . . . . .	36
Long ceiling installation cables . . . . .	19	Masimo rainbow SET single-patient use sensors . . . . .	36
Miscellaneous hardware . . . . .	20	Masimo rainbow SET reusable sensors . . . . .	36
Mounts for workstation . . . . .	20	<b>SpO2 monitoring – Masimo rainbow SET (RD rainbow SET)</b> . . . . .	37
Mounts for power supply . . . . .	20	Masimo rainbow RD-SET intermediate cables . . . . .	38
Mounts for the Cockpit and M500/M540 . . . . .	20	Masimo RD rainbow SET single-patient use sensors . . . . .	38
<b>IACS / M540 options</b> . . . . .	21	<b>SpO2 monitoring – Masimo rainbow SET (LNCS)</b> . . . . .	39
<b>ECG monitoring – reusable leadsets</b> . . . . .	22	Masimo rainbow SET LNCS intermediate cables . . . . .	40
<b>ECG monitoring – single-patient use leadsets</b> . . . . .	24	Masimo single-patient use LNCS sensors . . . . .	40
<b>ECG monitoring – MonoLead</b> . . . . .	26	Masimo reusable LNCS sensors . . . . .	41
ECG accessories . . . . .	26	Masimo single-patient use LNCS trauma sensors . . . . .	41
MonoLeads . . . . .	27	<b>SpO2 monitoring – Masimo rainbow SET (RD-SET)</b> . . . . .	42
<b>ECG monitoring – neonatal</b> . . . . .	28	Masimo rainbow SET MCable . . . . .	42
Neonatal leadsets and accessories . . . . .	28	Masimo rainbow RD-SET intermediate cables . . . . .	43
<b>Invasive pressure (IP) with Hemo2 Pod or Hemo4 pod</b> . . . . .	29	Masimo single-patient use RD-SET sensors . . . . .	43
Hemodynamic pods and cables . . . . .	29		
Transducer adapter cables . . . . .	29		
<b>IP Monitoring with MPod – QuadHemo or MCable – Dual Hemo</b> . . . . .	30		
IP cables . . . . .	30		
Transducer adapter cables . . . . .	30		
Cables for a second QuadHemo . . . . .	31		


Masimo single-patient use RD-SET Trauma sensors . . . . .	44	Temperature cables . . . . .	59
<b>SpO2 monitoring – Masimo SET (LNCS)</b> . . . . .	45	General purpose probes . . . . .	59
Masimo SET MCable . . . . .	45	Tcore accessories . . . . .	61
Masimo SET LNCS intermediate cables . . . . .	46	Tcore mounting accessories . . . . .	61
Masimo single-patient use LNCS sensors . . . . .	46	<b>MCable – Analog/Sync</b> . . . . .	62
Masimo reusable LNCS sensors . . . . .	47	Analog/sync cables . . . . .	62
Masimo single-patient use LNCS trauma sensors . . . . .	47	<b>Scio monitoring</b> . . . . .	63
<b>SpO2 monitoring – Masimo SET (RD SET)</b> . . . . .	48	Scio cables . . . . .	63
Masimo SET MCable . . . . .	48	Scio modules . . . . .	63
Masimo RD-SET intermediate cables . . . . .	49	Scio cables (M500 connection) . . . . .	64
Masimo single-patient use RD-SET sensors . . . . .	49	<b>NMT monitoring</b> . . . . .	65
Masimo reusable RD-SET sensors . . . . .	50	TOFscan module and sensor cables . . . . .	65
Masimo single-patient use RD-SET Trauma sensors . . . . .	50	TOFscan chargers . . . . .	65
<b>Masimo SpO2 accessories</b> . . . . .	51	TOFscan mounts . . . . .	65
<b>SpO2 monitoring – Nellcor Oximax</b> . . . . .	52		
Nellcor SpO2 cables . . . . .	52		
Nellcor single-patient use sensors . . . . .	52		
Nellcor reusable sensor . . . . .	53		
Nellcor FLEXMAX reusable sensors . . . . .	53		
<b>Cardiac (C.O.) output monitoring</b> . . . . .	54		
C.O. cables, catheters, and thermistors . . . . .	54		
Cardiac output kits (not shown) . . . . .	54		
<b>CO2 monitoring – Mainstream</b> . . . . .	55		
CO2 cables . . . . .	55		
Disposable airway adapters . . . . .	55		
Reusable airway adapters . . . . .	55		
Additional accessories (not shown) . . . . .	55		
<b>CO2 monitoring – Microstream</b> . . . . .	56		
CO2 cables . . . . .	56		
Sample lines for intubated patients . . . . .	56		
Sample lines for non-intubated patients . . . . .	57		
Mounting . . . . .	57		
<b>Temperature monitoring</b> . . . . .	58		
Temperature cables . . . . .	58		
Skin temperature sensors . . . . .	58		

## Infinity Acute Care System accessories overview

The following diagram is an overview of the accessories available with the Infinity Acute Care System.



Microstream CO2	Mainstream CO2	CO2 & Analog/Sync	Analog & Sync	Temp	
 <p>+</p> <p>+</p>	 <p>+</p> <p>Pediatric Adult</p>	 <p>+</p> <p>+</p> <p>CO2 Analog/Sync</p>	 <p>+</p> <p>+</p>	 <p>++</p> <p>+</p>	
ECG – reusable	ECG – MonoLead	ECG – single-use	ECG – neo	NIBP	
			 <p>+</p>	<p>Adult/Child Neonatal</p>  <p>+</p> <p>+</p>	



This chapter lists the accessories available with the Infinity Acute Care System (IACS).

The chapter is divided into the following sections:

- "General IACS accessories" on page 17
- "IACS / M540 options" on page 21
- "ECG monitoring – reusable leadsets" on page 22
- "ECG monitoring – single-patient use leadsets" on page 24
- "ECG monitoring – MonoLead" on page 26
- "ECG monitoring – neonatal" on page 28
- "Invasive pressure (IP) with Hemo2 Pod or Hemo4 pod" on page 29
- "IP Monitoring with MPod – QuadHemo or MCable – Dual Hemo" on page 30
- "Non-Invasive blood pressure (NIBP) monitoring – Dräger connectors" on page 32
- "Non-invasive blood pressure (NIBP) monitoring – neonatal" on page 34
- "SpO2 monitoring – Masimo rainbow SET (rainbow sensors)" on page 35
- "SpO2 monitoring – Masimo rainbow SET (RD rainbow SET)" on page 37
- "SpO2 monitoring – Masimo rainbow SET (LNCS)" on page 39
- "SpO2 monitoring – Masimo rainbow SET (RD-SET)" on page 42
- "SpO2 monitoring – Masimo SET (LNCS)" on page 45
- "SpO2 monitoring – Masimo SET (RD SET)" on page 48
- "Masimo SpO2 accessories" on page 51
- "SpO2 monitoring – Nellcor OxiMax" on page 52
- "Cardiac (C.O.) output monitoring" on page 54
- "CO2 monitoring – Mainstream" on page 55
- "CO2 monitoring – Microstream" on page 56

- "Temperature monitoring" on page 58
- "MCable – Analog/Sync" on page 62
- "Scio monitoring" on page 63
- "NMT monitoring" on page 65

On the following pages, each parameter section includes a diagram depicting parameter-specific accessories, along with one or more tables displaying accessory names and part numbers. The letters displayed in the diagrams correspond to the letters in the tables. Because some accessories are too numerous to list in the diagram, they do not have a letter associated with them.

Part numbers are subject to change without prior notice. Contact a Dräger representative for the most current part numbers.

Failure to use Dräger-approved accessories may compromise the functional integrity of the IACS.



## General IACS accessories

The following table lists the system components, power cables, system cables, and keyboards for the IACS.

### System components

The following parts are components in the IACS. They are not available for purchase individually. Contact a Dräger representative for additional information.

Part	Part No.
Infinity M500	MS20407
Infinity M540	MS20401
Infinity P2500	MS22277
Desktop 120W power supply (PS120)	2606270
IACS Y-piece for PS120	MS29702

### Cables and mouse

Part	Part No.
Mouse	MS23352
CMPACC Seal Shield mouse	MS27252
Mouse pad	MS24519
System cable, 1.5 m (4.9 ft)	MS20135
System cable, 2.5 m (8.2 ft)	MK08707
System cable, 3.5 m (11.5 ft)	MK08708
System cable, 4.5 m (14.7 ft)	MS20614
System cable, 2.2 m (7.2 ft) <i>(Note that this cable is only to be used in combination with Zeus IE)</i>	MK04408
Ethernet network cable, 4.9 m (16.1 ft)	4726399
Ethernet cables, 200.0 m (656.2 ft)	4725524
Ethernet cables, 165.0 m (541.3 ft)	4725532
Ethernet cables, 2.4 m (7.9 ft)	4726381
Ethernet cables, 1.2 m (3.9 ft)	4726373
Nurse call cable	8417370

Export protocol cable	MS22948
Patch cable, 20.0 m (65.6 ft)	4725557
Crossover cable, 6.0 m (19.7 ft)	5592667
Crossover cable, 2.0 m (6.6 ft)	5730176

### Recorders and printers

Part	Part No.
Infinity R50-N Recorder	5740068
Recorder Paper 50 mm	4711201
R50 Countertop Plate	5197384
Laser Printer 115V HP	6556513
Laserjet Printer HP 600, 220V	MS30021

### Power supplies

Part	Part No.
Power cable (North America)	4321720
Power cable (NEMA)	1868160
Power cable (Europe)	4321712
Power cable (Great Britain)	1851713
Power cable (Switzerland)	1851691
Power cable (Australia/New Zealand)	1851705
Power cable (China)	1859714
Power cable (Denmark)	1868950
Power cable (Israel)	1869833
Power cable (Brazil)	1875523
Power cable (Argentina)	1886274

## Keyboards

Part	Part No.
Washable keyboard (Russian)	MS24520
Washable keyboard (Croatian)	MS31000
Washable keyboard (Czech)	MS22286
Washable keyboard (Denmark)	MS20790
Washable keyboard (English)	MS22288
Washable keyboard (French)	MS22295
Washable keyboard (German)	MS22287
Washable keyboard (Greek)	MS23692
Washable keyboard (Hungarian)	MS22282
Washable keyboard (Italian)	MS22283
Washable keyboard (Norwegian)	MS22284
Washable keyboard (Portuguese)	MS22291
Washable keyboard (Romanian)	MS30998
Washable keyboard (Spanish)	MS22294
Washable keyboard (Swedish/ Finnish)	MS22285
Washable keyboard (Turkish)	MS22292

## Network panels and switches

Part	Part No.
Network patch panel	4725573
3Com network switch	7269850

## Device interface kits

The following device interface cables are available for the Infinity Acute Care System and are used to connect to the corresponding devices.

Device	Part No. (cable)
Edwards EV1000/Vigi.II/Vigileo	MS34114
Dräger Evita 2D/4/XL ventilator	
Dräger Savina 300 ventilator	
Dräger Carina ventilator	
Bis Vista monitor	MS34115
Dräger Primus/IE/Apollo	
Dräger Fabius Family	
Dräger V500 ventilator	MS34116
Dräger Babylog VN500	
Dräger V300 ventilator	
Dräger Oxylog 3000+ ventilator	
Dräger Perseus A500	
Dräger Zeus IE	
Maquet Servo-I ventilator	MS34117
TOFscan monitor	MS34118
TOF-Watch SX	N/A

Part	Part No.
Cockpit to DIM cable, 1.5 m (4.9 ft)	MS34112
Cockpit to DIM cable, 4.3 m (14.1 ft)	MS34113

## Long ceiling installation cables

Part	Part No.
Universal set, 2 x 5 m (6.5 x 16.4 ft)	G26455
Universal set, 2 x 8 m (6.5 x 26.25 ft)	G26456
IACS wall box assembly	G26740
IACS ceiling connection set	G25955

### Miscellaneous hardware

Part	Part No.
Transport dock and clamp	MS28144
Cable hook	MS26297
SpO <sub>2</sub> MCable mount	MS26266
Masimo rainbow SET MCable mount	MS28576
Rolling trolley for M540	RS-0025-05C

### Mounts for the Cockpit and M500/M540

Part	Part No.
Dräger lift arm, front rail, VESA, +300	G96548
Dräger double arm, front rail, VESA	G96715
M500 tilting mount	MS20427
Cockpit/M500 mount, VESA 75	MS20408
Mount spacer for M500	MS26141

### Mounts for workstation

Part	Part No.
Mount for workstation	MS17139
Mount bracket rail GCX	MS15934

#### NOTE

All other hardware from GCX can only be purchased directly from GCX corporation or one of its distributors.

### Mounts for power supply

Part	Part No.
Wall mounting adapter CPS/IDS power supply	4720061
Rail mount for CPS/IDS power supply	4720095
Wall bracket CPS/IDS power supply	4720129
Power supply rollstand bracket CPS/IDS	4720103
Mounting bracket for desktop power supply, 120 W	MS31216
Channel mount for power supply	MM15287

## IACS / M540 options

---


The following table lists the available IACS options. Check with a local sales representative for options that might be available for free during certain

ongoing promotions. The menu selections on the M540 or the C700/C500 are not available until the corresponding options have been enabled.

Option	M540 standalone	C700	C500	Part No.
M540 wireless	Yes	No	No	MS16266
Full arrhythmia	Yes	Yes	Yes	MS22225
Multi IP	Yes	Yes	Yes	MS20504
12-lead ECG	Yes	Yes	Yes	MS20508
12 channels	No	No	Yes	MS20505
16 channels	No	Yes	Yes	MS20506
Web-enabled layouts	No	Yes	Yes	MS20511
Web-enabled tabs	No	Yes	Yes	MS20512
View editor	No	Yes	Yes	MS20515
Physiological calculations	No	Yes	Yes	MS20516
Rainbow SpMet	Yes	Yes	Yes	MS28968
Rainbow SpCO	Yes	Yes	Yes	MS28969
Rainbow PVI	Yes	Yes	Yes	MS28970
Rainbow SpHb	Yes	Yes	Yes	MS28967
Rainbow SpOC	Yes	Yes	Yes	MS28972
Independent display*	No	Yes	Yes	MS33832

\* Not available on first-generation C500/C700 Cockpits or second-generation C500/C700 Cockpits

## ECG monitoring – reusable leadsets



### ECG accessories

	Part	Part No.
A	ECG extension cable 3/5/6, 2 m (6.5 ft) <sup>1)</sup>	MS16256
B	ECG extension cable, 4 lead, 2 m (6.5 ft)	MS31435
C	ECG 4-lead spacer, 10 pieces	MS22832

#### NOTE

<sup>1)</sup> The extension cable (MS16256) works with certain MonoLeads, traditional leads, and single-patient use leads. The extension cable is not compatible with ECG 4-lead grabber sets.


#### WARNING

To reduce the hazard of burns during electrosurgery, 4-lead extension cables are not supported for use in electrosurgical environments.

**Traditional reusable leadsets**

	<b>Part</b>	<b>Part No.</b>	<b>Old Part No.</b>
<b>D</b>	ECG grabber 3-lead, IEC1 (Europe), 1 m (3.3 ft)	MP03401	MS16159
	ECG grabber 3-lead, IEC2 (AHA/US), 1 m (3.3 ft)	MP03402	MS16231
<b>E</b>	ECG grabber 5-lead, IEC1 (Europe), 1 m (3.3 ft)	MP03403	MS16158
	ECG grabber 5-lead, IEC2 (AHA/US), 1 m (3.3 ft)	MP03404	MS16546
<b>F</b>	ECG grabber 6-lead, IEC1 (Europe), 1 m (3.3 ft)	MP03405	MS16157
	ECG grabber 6-lead, IEC2 (AHA/US), 1 m (3.3 ft)	MP03406	MS16547
<b>G</b>	ECG grabber 4-lead, IEC1 (Europe), 1 m (3.3 ft)	MP03407	MS22365
	ECG grabber 4-lead, IEC2 (AHA/US), 1 m (3.3 ft)	MP03408	MS22364

## ECG monitoring – single-patient use leadsets



### ECG accessories

	Part	Part No.
A	ECG extension cable 3/5/6, 2 m (6.5 ft) <sup>1)</sup>	MS16256
B	ECG extension cable, 4 lead, 2 m (6.5 ft)	MS31435
C	ECG 4-lead spacer, 10 pieces	MS22832

#### NOTE

1) The extension cable (MS16256) works with certain MonoLeads, traditional leads, and single-patient use leads. The extension cable is not compatible with ECG 4-lead grabber sets.

#### WARNING


To reduce the hazard of burns during electrosurgery, 4-lead extension cables are not supported for use in electrosurgical environments.



**Single-patient use leadsets**

	<b>Part</b>	<b>Part No.</b>
<b>D</b>	ECG grabber 3-lead, IEC1 (Europe), 1 m (3.3 ft)	MP00875
	ECG grabber 3-lead, IEC2 (AHA/US), 1 m (3.3 ft)	MP00877
<b>E</b>	ECG grabber 5-lead, IEC1 (Europe), 1 m / 1.5 m (3.3 ft / 4.9 ft)	MP00879
	ECG grabber 5-lead, IEC2 (AHA/US), 1 m / 1.5 m (3.3 ft / 4.9 ft)	MP00881
<b>F</b>	ECG grabber 6-lead, IEC1 (Europe), 1 m / 1.5 m (3.3 ft / 4.9 ft)	MP03122
	ECG grabber 6-lead, IEC2 (AHA/US), 1 m / 1.5 m (3.3 ft / 4.9 ft)	MP03123
<b>G</b>	ECG grabber 4-lead, IEC1 (Europe), 1 m (3.3 ft)	MP03120
	ECG grabber 4-lead, IEC2 (AHA/US), 1 m (3.3 ft)	MP03121

## ECG monitoring – MonoLead



### ECG accessories

	Part	Part No.
<b>A</b>	ECG extension cable 3/5/6, 2 m (6.5 ft) <sup>1)</sup>	MS16256
<b>B</b>	ECG extension cable, 4 lead, 2 m (6.5 ft) <sup>2)</sup>	MS31435
<b>C</b>	ECG 4-lead spacer, 10 pieces	MS22832

#### NOTE

<sup>1)</sup> The extension cable (MS16256) works with certain MonoLeads, traditional leads, and single-patient use leads. The extension cable is not compatible with ECG 4-lead grabber sets.

<sup>2)</sup> The 4-lead extension cable MS31435 is not for use with the 3.7 m (12.1 ft) long Monoleads MS20475 and MS20476.

#### WARNING

To reduce the hazard of burns during electrosurgery, 4-lead extension cables are not supported for use in electrosurgical environments.

## MonoLeads

	Part	Part No.
<b>D</b>	ECG MonoLead 3-lead, IEC1 (Europe), 2 m (6.5 ft)	MS14555
	ECG MonoLead 3-lead, IEC2 (AHA/US), 2 m (6.5 ft)	MS14556
	ECG MonoLead 3-lead, IEC1 (Europe), 4.1 m (13.5 ft) <sup>1)</sup>	MS28561
	ECG MonoLead 3-lead, IEC2 (AHA/US), 4.1 m (13.5 ft) <sup>1)</sup>	MS28557
<b>E</b>	ECG MonoLead 5-lead, IEC1 (Europe), 2.5 m (8.2 ft)	MS14559
	ECG MonoLead 5-lead, IEC2 (AHA/US) 2.5 m (8.2 ft)	MS14560
	ECG MonoLead 5-lead, IEC1 (Europe), 5.3 m (17.4 ft) <sup>1)</sup>	MS28559
	ECG MonoLead 5-lead, IEC2 (AHA/US), 5.3 m (17.4 ft) <sup>1)</sup>	MS28558
<b>F</b>	ECG MonoLead 6-lead, IEC1 (Europe), 3 m (10 ft)	MS14683
	ECG MonoLead 6-lead, IEC2 (AHA/US), 3 m (10 ft)	MS14682
	ECG MonoLead 6-lead, IEC1 (Europe), 5.7 m (18.7 ft) <sup>1)</sup>	MS17185
	ECG MonoLead 6-lead, IEC2 (AHA/US), 5.7 m (18.7 ft) <sup>1)</sup>	MS28560
<b>G</b>	ECG MonoLead 4-lead, IEC1 (Europe), 1.7 m (5.6 ft)	MS14693
	ECG MonoLead 4-lead, IEC2 (AHA/US), 1.7 m (5.6 ft)	MS14694
	ECG MonoLead 4-lead, IEC1 (Europe), 3.7 m (12.1 ft) <sup>2)</sup>	MS20476
	ECG MonoLead 4-lead, IEC2 (AHA/US), 3.7 m (12.1 ft) <sup>2)</sup>	MS20475

### NOTE

<sup>1)</sup>The indicated accessories must *not* be used in combination with the extension cable MS16256.


<sup>2)</sup>The indicated accessories must *not* be used in combination with the extension cable MS31435.

## ECG monitoring – neonatal

### Neonatal leadsets and accessories


#### WARNING

To prevent burns during electrosurgery (ESU), do not use the ECG neonatal 3-lead adapter cable with the ECG extension cable (MS16256).



	Part	Part No.
<b>A</b>	ECG neonatal 3-lead adapter cable, 2 m (6.5 ft), single-pin	MS25951
<b>B</b>	ECG electrode KittyKat 300 pieces	5195024
<b>C</b>	ECG 4-lead spacer, 10 pieces	MS22832

## Invasive pressure (IP) with Hemo2 Pod or Hemo4 pod



### Hemodynamic pods and cables

	Part	Part No.
A	M540 Hemo Connect cable, 1.5 m (4.9 ft) <sup>1)</sup>	MS23897
	M540 Hemo Connect cable, 2.9 m (9.5 ft) <sup>1)</sup>	MS23898
B	Hemo2 pod ( <i>not shown</i> – no longer available)	4319435
	Hemo4 pod (no longer available)	4315961


#### NOTE

<sup>1)</sup> The M540 Hemo Connect cable is used with the Hemo2 pod/Hemo4 pod.

### Transducer adapter cables

	Part	Part No.
C	Hemo pod adapter, 10-pin	3375958
D	Hemo pod adapter, SensoNor, 7-pin	MS23598
E	Hemo pod adapter, Baxter-Edwards	5196980
F	Hemo pod adapter, Abbott/Medex/ICU Medical <sup>®</sup>	5196998
G	Hemo pod adapter, Argon/Merit Medical <sup>®</sup>	3375941

## IP Monitoring with MPod – QuadHemo or MCable – Dual Hemo



### IP cables

	Part	Part No.
<b>A</b>	MCable – Dual Hemo	MS20783
<b>B</b>	M540 - Hemo Connect cable, 1.5 m (4.9 ft) <sup>1)</sup>	MS23897
	M540 Hemo Connect cable, 2.9 m (9.5 ft) <sup>1)</sup>	MS23898
<b>C</b>	MPod – QuadHemo	MS20725


#### NOTE

<sup>1)</sup> The M540 Hemo Connect cable is used with the MPod – QuadHemo (not necessary for the MCable – Dual Hemo).

### Transducer adapter cables

	Part	Part No.
<b>D</b>	Abbott/Medex/ICU Medical®	MS22535
<b>E</b>	Baxter-Edwards	MS22147
<b>F</b>	Argon/Merit Medical®	MS22148
<b>G</b>	Utah Medical	MS22534
<b>H</b>	Dräger 7-pin	MS22533
<b>I</b>	Dräger 10-pin	MS22532

### Cables for a second QuadHemo




	Part	Part No.
A	Pod communication cable, 0.3 m (1.0 ft)	7257988
	Pod communication cable, 3.0 m (9.8 ft)	3368425
	Pod communication cable, 5.0 m (16.4 ft)	5195198

#### NOTE

The pod communication cable is used to daisy chain two MPod – QuadHemo devices for simultaneous use.

## Non-Invasive blood pressure (NIBP) monitoring – Dräger connectors



### NIBP hose

The following single-patient use and reusable cuffs and hose use Dräger connectors.

	Part	Part No.
<b>A</b>	Connection hose, adult/pediatric 3.7 m (12.1 ft)	MP00953

### Reusable adult/pediatric blood pressure cuffs


#### NOTE

All cuffs are not made with natural rubber latex.

	Part	Part No.
<b>B</b>	Infant arm cuff, XXS 8 to 13 cm (3 to 5 in)	MP00911
	Child arm cuff, XS 12 to 19 cm (5 to 7 in)	MP00912
	Small adult arm cuff, S 17 to 29 cm (6 to 10 in)	MP00913
	Adult arm cuff, M 23 to 33 cm (9 to 13 in)	MP00915
	Adult arm cuff, M+ 23 to 43 cm (9 to 17 in)	MP00916
	Adult arm cuff, M++ 23 to 53 cm (9 to 21 in)	MP00917
	Large adult arm cuff, L 31 to 40 cm (12 to 16 in)	MP00918
	Large adult arm cuff, L+ 31 to 55 cm (12 to 22 in)	MP00919
	Large adult thigh cuff, XL 38 to 50 cm (15 to 20 in)	MP00921




**Single-patient use adult/pediatric blood pressure cuffs**



	<b>Part</b>	<b>Part No.</b>
<b>B</b>	Infant arm cuff, XXS, 8 to 13 cm (3 to 5 in) 10 pieces	MP00924
	Child arm cuff, XS, 12 to 19 cm (5 to 7 in) 10 pieces	MP00925
	Small adult arm cuff, S, 17 to 29 cm (6 to 11 in) 10 pieces	MP00926
	Adult arm cuff, M, 23 to 33 cm (9 to 13 in) 10 pieces	MP00928
	Adult arm cuff, M+, 23 to 43 cm (9 to 17 in) 10 pieces	MP00929
	Adult arm cuff, M++, 23 to 53 cm (9 to 21 in) 10 pieces	MP00930
	Adult arm cuff, L, 31 to 40 cm (12 to 16 in) 10 pieces	MP00931
	Adult arm cuff, L+, 31 to 55 cm (12 to 22 in) 10 pieces	MP00932
	Adult thigh cuff, XL, 38 to 50 cm (15 to 20 in) 10 pieces	MP00934

## Non-invasive blood pressure (NIBP) monitoring – neonatal



### NIBP hose


The following single-patient use and reusable cuffs and hoses use universal connectors.

	Part	Part No.
<b>A</b>	Neonatal NIBP hose, 2.4 m (7.87 ft)	2870298

### Single-patient use neonatal blood pressure cuffs

	Part (arm circumference)	Part No.
<b>B</b>	Neonatal cuff #1 3.1 to 5.7 cm (1.2 to 2.2 in)	2870181
	Neonatal cuff #2 4.3 to 8.0 cm (1.7 to 3.1 in)	2870199
	Neonatal cuff #3 5.8 to 10.9 cm (2.3 to 4.3 in)	2870207
	Neonatal cuff #4 7.1 to 13.1 cm (2.8 to 5.1 in)	2870215
	Neonatal cuff #5 8.3 to 15.0 cm (3.2 to 6 in)	2870173

## SpO2 monitoring – Masimo rainbow SET (rainbow sensors)



### Masimo rainbow SET MCable

	Part	Part No.
A	Infinity MCable – Masimo rainbow SET	MS27003*
* Order using part number MS27900		

### Masimo rainbow SET intermediate cables

	Part	Part No. (non-US)	Part No. (US)	Part No. (Japan)
<b>B</b>	Masimo rainbow cable RC-4, 1.2 m (4 ft)	MP02974	MP03052	MP03154
	Masimo rainbow cable RC-12, 3.6 m (12 ft)	MP02975	MP03053	MP03155

### Masimo rainbow SET single-patient use sensors

Measures SpCO, SpMet, SpO<sub>2</sub>, SpOC, PI, PLS,  
PVI <sup>1)</sup>.

	Part	Part No. (non-US)	Part No. (US)	Part No. (Japan)
<b>C</b>	Masimo rainbow sensor R25-L, 10 pieces	MP02961	MP03054	MP03156
	Masimo rainbow sensor R20-L, 10 pieces	MP02962	MP03055	MP03157
	Masimo rainbow sensor R25, 10 pieces	MP02959	MP03056	MP03158
	Masimo rainbow sensor R20, 10 pieces	MP02960	MP03057	MP03159

**NOTE**

<sup>1)</sup> PVI works only with Masimo rainbow SET.

### Masimo rainbow SET reusable sensors


Measures SpCO, SpMet, SpO<sub>2</sub>, SpOC, PI, PLS,  
PVI <sup>1)</sup>.

	Part	Part No. (non-US)	Part No. (US)	Part No. (Japan)
<b>D</b>	Masimo rainbow sensor DCI, adult	MP02971	MP03064	MP03166
	Masimo rainbow sensor DCIP, pediatric	MP02972	MP03065	MP03167

**NOTE**

<sup>1)</sup> PVI works only with Masimo rainbow SET.

## SpO2 monitoring – Masimo rainbow SET (RD rainbow SET)



### Masimo rainbow SET MCable

	Part	Part No.
<b>A</b>	Infinity MCable – Masimo rainbow SET	MS27003*
* Order using part number MS27900		

### Masimo rainbow RD-SET intermediate cables

	Part	Part No. (non-US)	Part No. (US)	Part No. (Japan)
<b>B</b>	Masimo Cbl RD rainbow SET 1.5m (4.9 ft)	MS33738	MS33757	MS33776
	SpO2 Masimo Cbl RD rainbow SET 3.6m (11.8 ft)	MS33739	MS33758	MS33777

**NOTE**

<sup>1)</sup> PVI works only with Masimo rainbow SET.

### Masimo RD rainbow SET single-patient use sensors


Measures SpHb, SpMet, SpO<sub>2</sub>, SpOC, PI, PLS,  
PVI <sup>1)</sup>.

	Part	Part No. (non-US)	Part No. (US)	Part No. (Japan)
<b>C</b>	SpO2 Masimo RD rainbow SET, adult, 10 pieces	MS40076	MS40074	MS40075
	SpO2 Masimo RD rainbow SET, pediatric, 10 pieces	MS40079	MS40077	MS40078
<b>D</b>	SpO2 Masimo RD rainbow SET, infant, 10 pieces	MS40082	MS40080	MS40081
	SpO2 Masimo RD rainbow SET, neonatal, 10 pieces	MS40085	MS40083	MS40084

**NOTE**

<sup>1)</sup> PVI works only with Masimo rainbow SET.

## SpO2 monitoring – Masimo rainbow SET (LNCS)



### Masimo rainbow SET MCable

	Part	Part No.
A	Infinity MCable – Masimo rainbow SET	MS27003*
* Order using part number MS27900		

**Masimo rainbow SET LNCS intermediate cables**

	<b>Part</b>	<b>Part No. (non-US)</b>	<b>Part No. (US)</b>	<b>Part No. (Japan)</b>
<b>B</b>	SpO2 Masimo red LNC-4 cable, 1.2 m (4 ft)	MP02991	MP03030	MP03109
	SpO2 Masimo red LNC-10 cable, 3 m (10 ft)	MP02992	MP03031	MP03110
	SpO2 Masimo red LNC-14 cable, 4.3 m (14 ft)	MP03196	MP03032	MP03111

**Masimo single-patient use LNCS sensors**

	<b>Part</b>	<b>Part No. (non-US)</b>	<b>Part No. (US)</b>	<b>Part No. (Japan)</b>
<b>C</b>	SpO2 Masimo LNCS, adult, 20 pcs	MP00790	MP03011	MP03089
	SpO2 Masimo LNCS, infant, 20 pcs	MP00791	MP03013	MP03091
	SpO2 Masimo LNCS, neonatal, 20 pcs	MP00792	MP03014	MP03092
	SpO2 Masimo LNCS, preterm neonatal, 20 pcs	MP00794	MP03015	MP03093
	SpO2 Masimo LNCS, pediatric, 20 pcs	MP00793	MP03012	MP03090



**Masimo reusable LNCS sensors**

	<b>Part</b>	<b>Part No. (non-US)</b>	<b>Part No. (US)</b>	<b>Part No. (Japan)</b>
<b>D</b>	SpO2 Masimo LNCS, DCI, adult	MP00796	MP03019	MP03097
	SpO2 Masimo, LNCS, DCIP, pediatric	MP00795	MP03020	MP03098
	TC-I adult ear	MP00788	MP03021	MP03099
	TF-I adult forehead	MP00799	MP03022	MP03101
	YI adult/pediatric/ infant/neonatal multi-site	MP00789	MP03023	MP03102
	Masimo LNCS-DBI, adult <sup>1)</sup>	MP02993	MP03024	MP03103

**NOTE**

<sup>1)</sup> Masimo does not make any claims regarding motion.

**Masimo single-patient use LNCS trauma sensors**

	<b>Part</b>	<b>Part No. (non-US)</b>	<b>Part No. (US)</b>	<b>Part No. (Japan)</b>
<b>E</b>	Masimo LNCS trauma sensor, adult	MP03191	MP03016	MP03094
	Masimo LNCS trauma sensor, neonatal	MP03192	MP03017	MP03095
	Masimo LNCS trauma sensor, infant/pediatric	MP03193	MP03018	MP03096


**CAUTION**

With trauma sensors, the sensor-off detection performance may be compromised.

## SpO2 monitoring – Masimo rainbow SET (RD-SET)

### Masimo rainbow SET MCable

	Part	Part No.
A	Infinity MCable – Masimo rainbow SET	MS27003*
* Order using part number MS27900		



576

**Masimo rainbow RD-SET intermediate cables**

	Part	Part No. (non-US)	Part No. (US)	Part No. (Japan)
<b>B</b>	Masimo Cbl RD rainbow SET 1.5m (4.9 ft)	MS33738	MS33757	MS33776
	SpO <sub>2</sub> Masimo Cbl RD rainbow SET 3.6m (11.8 ft)	MS33739	MS33758	MS33777

**Masimo single-patient use RD-SET sensors**

	Part	Part No. (non-US)	Part No. (US)	Part No. (Japan)
<b>C</b>	SpO <sub>2</sub> Masimo RD-SET adult, 20 pieces	MS33726	MS33745	MS33764
	SpO <sub>2</sub> Masimo RD-SET pediatric, 20 pieces	MS33727	MS33746	MS33765
	SpO <sub>2</sub> Masimo RD-SET infant, 20 pieces	MS33728	MS33747	MS33766
	SpO <sub>2</sub> Masimo RD-SET neonatal, 20 pieces	MS33729	MS33748	MS33767
	SpO <sub>2</sub> Masimo RD-SET preterm, 20 pieces	MS33730	MS33749	MS33768

**Masimo reusable RD-SET sensors**

	Part	Part No. (non-US)	Part No. (US)	Part No. (Japan)
<b>D</b>	SpO <sub>2</sub> Masimo RD-SET DCI, adult	MS33740	MS33759	MS33778
	SpO <sub>2</sub> Masimo RD-SET DCIP, pediatric	MS33741	MS33760	MS33779
	SpO <sub>2</sub> Masimo RD-SET TC-I, adult ear	MS40067	MS40065	MS40066
	SpO <sub>2</sub> Masimo RD-SET TF-I, adult forehead	MS40073	MS40071	MS40072
	SpO <sub>2</sub> Masimo RD-SET YI, adult/pediatric/infant/neonatal multi-site	MS40070	MS40068	MS40069


**Masimo single-patient use RD-SET Trauma sensors**

	<b>Part</b>	<b>Part No. (non-US)</b>	<b>Part No. (US)</b>	<b>Part No. (Japan)</b>
<b>E</b>	SpO2 Masimo RD-SET trauma sensor, adult, 20 pieces	MS40088	MS40086	MS40087
	SpO2 Masimo RD-SET trauma sensor, infant/pediatric, 20 pieces	MS40091	MS40089	MS40090
	SpO2 Masimo RD-SET trauma sensor, neonatal, 20 pieces	MS40094	MS40092	MS40093

**CAUTION**

With trauma sensors, the sensor-off detection performance may be compromised.

## SpO2 monitoring – Masimo SET (LNCS)



### Masimo SET MCable

	Part	Part No.
A	Infinity MCable – Masimo SET	MS20667

**Masimo SET LNCS intermediate cables**

	<b>Part</b>	<b>Part No. (non-US)</b>	<b>Part No. (US)</b>	<b>Part No. (Japan)</b>
<b>B</b>	Masimo intermediate cable, LNCS 1.2 m (4 ft)	MP00800	MP03027	MP03106
	Masimo intermediate cable, LNCS 3 m (10 ft)	MP00797	MP03025	MP03107
	Masimo intermediate cable, LNCS 4.3 m (14 ft)	MP03194	MP03026	MP03108

**Masimo single-patient use LNCS sensors**

	<b>Part</b>	<b>Part No. (non-US)</b>	<b>Part No. (US)</b>	<b>Part No. (Japan)</b>
<b>C</b>	SpO2 Masimo LNCS, adult, 20 pcs	MP00790	MP03011	MP03089
	SpO2 Masimo LNCS, infant, 20 pcs	MP00791	MP03013	MP03091
	SpO2 Masimo LNCS, neonatal, 20 pcs	MP00792	MP03014	MP03092
	SpO2 Masimo LNCS, preterm neonatal, 20 pcs	MP00794	MP03015	MP03093
	SpO2 Masimo LNCS, pediatric, 20 pcs	MP00793	MP03012	MP03090

**Masimo reusable LNCS sensors**

	<b>Part</b>	<b>Part No. (non-US)</b>	<b>Part No. (US)</b>	<b>Part No. (Japan)</b>
<b>D</b>	SpO2 Masimo LNCS, DCI, adult	MP00796	MP03019	MP03097
	SpO2 Masimo LNCS, DCIP, pediatric	MP00795	MP03020	MP03098
	TC-I adult ear	MP00788	MP03021	MP03099
	TF-I adult forehead	MP00799	MP03022	MP03101
	YI adult/pediatric/ infant/neonatal multi-site	MP00789	MP03023	MP03102
	Masimo LNCS-DBI, adult <sup>1)</sup>	MP02993	MP03024	MP03103

**NOTE**

<sup>1)</sup> Masimo does not make any claims regarding motion.

**Masimo single-patient use LNCS trauma sensors**

	<b>Part</b>	<b>Part No. (non-US)</b>	<b>Part No. (US)</b>	<b>Part No. (Japan)</b>
<b>E</b>	Masimo LNCS trauma sensor, adult	MP03191	MP03016	MP03094
	Masimo LNCS trauma sensor, neonatal	MP03192	MP03017	MP03095
	Masimo LNCS trauma sensor, infant/pediatric	MP03193	MP03018	MP03096


**CAUTION**

With trauma sensors, the sensor-off detection performance may be compromised.

## SpO2 monitoring – Masimo SET (RD SET)

### Masimo SET MCable

	Part	Part No.
A	Infinity MCable – Masimo SET	MS20667



574



**Masimo RD-SET intermediate cables**

	Part	Part No. (non-US)	Part No. (US)	Part No. (Japan)
<b>B</b>	Masimo cable RD SET, 1.5m (4.9 ft)	MS33736	MS33755	MS33774
	Masimo cable RD SET, 3.6m (11.8 ft)	MS33737	MS33756	MS33775

**Masimo single-patient use RD-SET sensors**

	Part	Part No. (non-US)	Part No. (US)	Part No. (Japan)
<b>C</b>	SpO2 Masimo RD-SET adult, 20 pieces	MS33726	MS33745	MS33764
	SpO2 Masimo RD-SET pediatric, 20 pieces	MS33727	MS33746	MS33765
	SpO2 Masimo RD-SET infant, 20 pieces	MS33728	MS33747	MS33766
	SpO2 Masimo RD-SET neonatal, 20 pieces	MS33729	MS33748	MS33767
	SpO2 Masimo RD-SET preterm, 20 pieces	MS33730	MS33749	MS33768

**Masimo reusable RD-SET sensors**

	<b>Part</b>	<b>Part No. (non-US)</b>	<b>Part No. (US)</b>	<b>Part No. (Japan)</b>
<b>D</b>	SpO2 Masimo RD-SET DCI, adult	MS33740	MS33759	MS33778
	SpO2 Masimo RD-SET DCIP, pediatric	MS33741	MS33760	MS33779
	SpO2 Masimo RD-SET TC-I, adult ear	MS40067	MS40065	MS40066
	SpO2 Masimo RD-SET TF-I, adult forehead	MS40073	MS40071	MS40072
	SpO2 Masimo RD-SET YI, adult/pediatric/infant/neonatal multi-site	MS40070	MS40068	MS40069

**Masimo single-patient use RD-SET Trauma sensors**

	<b>Part</b>	<b>Part No. (non-US)</b>	<b>Part No. (US)</b>	<b>Part No. (Japan)</b>
<b>E</b>	SpO2 Masimo RD-SET trauma sensor, adult, 20 pieces	MS40088	MS40086	MS40087
	SpO2 Masimo RD-SET trauma sensor, infant/pediatric, 20 pieces	MS40091	MS40089	MS40090
	SpO2 Masimo RD-SET trauma sensor, neonatal, 20 pieces	MS40094	MS40092	MS40093

**CAUTION**


With trauma sensors, the sensor-off detection performance may be compromised.

## Masimo SpO<sub>2</sub> accessories

---

<b>Part</b>	<b>Part No. (non-US)</b>	<b>Part No. (US)</b>	<b>Part No. (Japan)</b>
Adhesive SpO <sub>2</sub> square, 12 pieces	7496875	MP03048	MP03150
Petite wrap-adult/pedi/neo	7496909	MP03047	MP03149
Posey wrap, 12 pieces	7496933	MP03045	MP03147
SpO <sub>2</sub> wrap for Masimo YI	7496917	MP03046	MP03148
SpO <sub>2</sub> foam wrap for Masimo YI (12 pieces)	7496891	MP03049	MP03151
Headband for LNOP/LNCS TF-I	MP03197	MP03050	MP03152
Adhesive pad LNOP/LNCS TF-I	MP03198	MP03051	MP03153
Masimo rainbow sensor shield	MP03203	MP03070	MP03172

## SpO<sub>2</sub> monitoring – Nellcor OxiMax



### NOTE

Purchase of this instrument confers no express or implied license under any Nellcor patent to use this instrument with any oximetry sensor that is not manufactured or licensed by Nellcor.

### Nellcor SpO<sub>2</sub> cables

	Part	Part No.
<b>A</b>	Infinity MCable – Nellcor OxiMax	MS20668
<b>B</b>	Nellcor intermediate cable, 1.2 m (3.9 ft)	MP00748
	Nellcor intermediate cable, 3 m (9.8 ft)	MS20979

### Nellcor single-patient use sensors

The following tables list the compatible sensors available for purchase from Dräger.

	Part	Part No.
<b>C</b>	OxiMax MaxA, 46 cm (18 in), 24 pieces	MX50065
	OxiMax MaxAL, 91 cm (36 in), 24 pieces	MX50071
	OxiMax MaxI, infant, 24 pieces	MX50067
	OxiMax MaxN, neo/adult, 24 pieces	MX50068
	OxiMax MaxP, pediatric, 24 pieces	MX50066
	Nellcor, OxiMax Max-R	MX50069
	Nellcor SpO <sub>2</sub> sensor, OxiMax Max-Fast, 24 pieces	MX50070
	OxiBand Oxi-P/I, ped/inf, 100 pieces	2170221
	Nellcor Oxiband wrap A/N, 100 pieces	8201012

**Nellcor reusable sensor**


	<b>Part</b>	<b>Part No.</b>
<b>D</b>	Durasensor DS100A, adult, finger	7262764
	Nellcor DURA-Y SpO2 sensor, D-YS	MX01004

**Nellcor FLEXMAX reusable sensors**

	<b>Part</b>	<b>Part No.</b>
<b>E</b>	Nellcor FLEXMAX Adt, SpO2 sensor, adult, finger	MS34677
	Nellcor FLEXMAX Pdt, SpO2 sensor, pediatric, finger	MS34678

## Cardiac (C.O.) output monitoring

### C.O. cables, catheters, and thermistors



	Part	Part No.
A	M540 Hemo Connect cable, 1.5 m (4.9 ft)	MS23897
	M540 Hemo Connect cable, 2.9 m (9.5 ft)	MS23898
B	MPod – QuadHemo	MS20725
	Hemo2 pod (not shown)	4319435
	Hemo4 pod (not shown)	4315961
C	Intermediate cable, 1 m (3.3 ft)	3368458
D	Intermediate cable to catheter	8419160
E	Thermistor cable, Argon/Merit Medical®	8420077
F	Thermistor T-piece, 25 pieces	5741975
G	Thermistor cable, Baxter-Edwards	8539983


### Cardiac output kits (not shown)

Part	Part No.
Cardiac output kit, Baxter <sup>1)</sup>	7265155
Cardiac output kit, Argon/Merit Medical® <sup>1)</sup>	7265163

#### NOTE

<sup>1)</sup> Kit includes intermediate cable, catheter cable, and thermistor cable.

## CO<sub>2</sub> monitoring – Mainstream



### CO<sub>2</sub> cables

	Part	Part No.
A	MCable - Mainstream CO <sub>2</sub>	6871950
B	Y-cable CO <sub>2</sub> and analog/sync <sup>1)</sup>	MS22259

#### NOTE

<sup>1)</sup> The Y-cable allows simultaneous use of the MCable - Mainstream CO<sub>2</sub> and MCable – Analog/Sync.

### Disposable airway adapters

	Part	Part No.
C	Airway adapter main, adult	MP01062
	Airway adapter main, pediatric	MP01063


### Reusable airway adapters

	Part	Part No.
C	Airway adapter main, adult	6870279
	Airway adapter main, pediatric	6870280

### Additional accessories (not shown)

Part	Part No.
Calibration set	8412710
CO <sub>2</sub> test gas cylinder	6850435
Parking holder for CO <sub>2</sub> sensor	8416352
Test filter	6870281

## CO2 monitoring – Microstream



### CO2 cables

	Part	Part No.
<b>A</b>	Infinity MCable – Microstream CO2	MS32641
<b>B</b>	Y-cable CO2 and analog/sync <sup>1)</sup>	MS22259
	MCable – Microstream CO2 extension cable	6872159

#### NOTE

<sup>1)</sup> The Y-cable allows simultaneous use of the MCable – Microstream and MCable – Analog/Sync.

### Sample lines for intubated patients

	Part	Part No.
<b>C</b>	etCO2 Airway adapter with sample line, adult/pediatric, max. 12 h, 25 pcs	7869535
	etCO2 Airway adapter with sample line, neonatal/pediatric, max. 72 h, 25 pcs	7869550
	etCO2 Airway adapter with sample line, adult/pediatric, max. 72 h, 25 pcs	7869543




## Sample lines for non-intubated patients

	Part	Part No.
<b>C</b>	etCO <sub>2</sub> nasal cannula, child, max. 12 h, 25 pcs	7869469
	etCO <sub>2</sub> nasal cannula, adult, max. 24 h, 25 pcs	7869592
	etCO <sub>2</sub> nasal cannula, + O <sub>2</sub> , child, max. 24 h, 25 pcs	7869485
	etCO <sub>2</sub> nasal cannula, + O <sub>2</sub> , adult, max. 24 h, 25 pcs	7869493
	etCO <sub>2</sub> nasal/oral cannula, + O <sub>2</sub> , child, max. 12 h, 25 pcs	7869584
	Smart CapnoLine Plus adult/intermediate, 25 pcs	MS26187
	Smart CapnoLine Plus O <sub>2</sub> adult/intermediate, 25 pcs	MS26188
	etCO <sub>2</sub> nasal cannula, neonate/child, max. 24 h, 25 pcs	7869618
	etCO <sub>2</sub> nasal/oral cannula, child, max. 12 h, 25 pcs	7869501
	etCO <sub>2</sub> nasal cannula, adult, max. 12 h, 25 pcs	7869477

## Mounting

Part	Part No.
Universal pole mount	7485621
MCable pole mount adapter	MS32774
MCable holder	MS32642

## Temperature monitoring



### Temperature cables

	Part	Part No.
A	Adapter cable 7-pin connector for temperature sensor with 1/4-inch connector	5198333
	Dual temperature Y cable, 7-pin	5592154

### Skin temperature sensors

#### WARNING

Skin temperature sensors are not for use in electrosurgery.

#### Adult reusable


	Part	Part No.
B	Skin temperature sensor, 1.5 m (4.9 ft)	4329822
	Skin temperature sensor, 3 m (9.8 ft)	5204669

#### Adult single-patient use

	Part	Part No.
C	Skin temperature sensor, adult, 1.6 m (5.2 ft)	MP00995
	Skin temperature sensor, adult, 3 m (9.8 ft)	MP00996

#### Pediatric/neonatal single-patient use

	Part	Part No.
C	Skin temperature sensor, pediatric, 1.6 m (5.2 ft)	MP00997
	Skin temperature sensor, pediatric, 3 m (9.8 ft)	MP00998



### Temperature cables

	Part	Part No.
A	Adapter cable 7-pin connector for temperature sensor with 1/4-inch connector	5198333
	Dual temperature Y cable, 7-pin	5592154

### General purpose probes

#### WARNING


The temperature sensor cover 7014616 contains latex.

#### NOTE

Cover internally placed reusable temperature sensors with temperature sensor sheaths.

### Adult reusable general purpose probes

	Part	Part No.
B	Temperature sensor, adult, 1.5 m (4.9 ft)	4329889
	Temperature sensor, adult, 3 m (9.8 ft)	5204644
C	Core probe covers, 20 cm, 10 pieces (contains latex)	7014616
	Temperature sensor, adult, 1.5 m (4.9 ft) and core probe covers, 20 cm, 10 pieces (contains latex)	7498939
	Core probe covers, 20 cm, 100 pieces (latex-free)	MS40243



**Pediatric reusable general purpose probes**

	<b>Part</b>	<b>Part No.</b>
<b>B</b>	Temperature sensor, pediatric, 1.5 m (4.9 ft)	4329848
	Temperature sensor, pediatric, 3 m (9.8 ft)	5204651
	Temperature sensor, pediatric, 1.5 m (4.9 ft) and core probe cover, 20 cm, 10 pieces (contains latex)	7498921

**Adult single-patient use general purpose probes**

<b>D</b>	Temperature sensor, adult, 1.6 m (5.2 ft)	MP00991
	Temperature sensor, adult, 3 m (9.8 ft)	MP00992

**Pediatric single-patient use general purpose probes**


<b>D</b>	Temperature sensor, pediatric, 1.6 m (5.2 ft)	MP00993
	Temperature sensor, pediatric, 3 m (9.8 ft)	MP00994

**NOTE**

<sup>1)</sup> Temperature core probes can be used during electrosurgery and in esophageal and rectal applications.

**NOTE**

Temperature sensors with 7-pin connectors are compatible with Hemo4 pod, Hemo2 pod, MPod – QuadHemo (see page 29).



**Tcore accessories**

	<b>Part</b>	<b>Part No.</b>
<b>A</b>	Tcore adapter (reusable)	MP00999
<b>B</b>	Tcore sensor (single-patient use), 20 pieces	MP00989
	Tcore adapter extension cable, 1.6 m (5.2 ft)	MP02201

**Tcore mounting accessories**

<b>Part</b>	<b>Part No.</b>
Holder for MultiMed	MP00721
Tcore mount for M540*	MS33239

\* Mounts on top of SpO2 MCable holder


## MCable – Analog/Sync

### Analog/sync cables

	Part	Part No.
<b>A</b>	Y-cable CO <sub>2</sub> and analog/sync <sup>1)</sup>	MS22259
<b>B</b>	MCable – Analog/Sync	MS20662
<b>C</b>	Analog output cable, 5 m (16.4 ft)	4314618
	QRS sync output cable, 3 m (9.8 ft)	4314667


#### NOTE

<sup>1)</sup> The y-cable allows simultaneous use of the MCable - Mainstream CO<sub>2</sub> and MCable – Analog/Sync.



127

## Scio monitoring



### Scio cables

	Part	Part No.
A	Scio connection cable 2.2 m (7.2 ft)	MS33167
	Scio connection cable 0.8 m (2.6 ft)	MS33168

### Scio modules


	Part	Part No.
B	Scio Four	6871804
	Scio Four plus	6871802
	Scio Four Oxi	6871803
	Scio Four Oxi plus	6871801

#### NOTE

See the Scio Four Modules supplement for additional accessories.

**Scio cables (M500 connection)**

	Part	Part No.
<b>C</b>	Scio to M500 connection cable 0.8 m (2.6 ft)	MS33986
	Scio to M500 connection cable 2.2 m (7.2 ft)	MS33987



578



## NMT monitoring

---

### TOFscan module and sensor cables

Part	Part No.
TOFscan device	MS32146*
Hand sensor and stimulation cable	MS32154
Pediatric sensor	MS32289
Eyebrow sensor	MS32291
Foot sensor	MS32290
Extension cable 1.8 m (5.9 ft)	MS32153
Optical adapter	MS32156
* Order using part number MS32318	

### TOFscan chargers

Part	Part No.
AC power / charger (EU)	MS32147
AC power / charger (UK)	MS32148
AC power / charger (US)	MS32149
AC power / charger (Type 1)	MS34186
AC power / charger (Universal, requires locally purchased cable for connection to power socket)	MS34187

### TOFscan mounts

Part	Part No.
Clamp fixation (diameter 10-40 mm)	MS32150
Clamp fixation (diameter 20-60 mm)	MS32151
Cable support	MS33884

This page has been left blank intentionally.






These instructions for use only apply to  
**Infinity® Acute Care System –  
Monitoring accessories VG7.n**  
with the Serial No.:



If no Serial No. has been filled in by Dräger,  
these Instructions for Use are provided for  
general information only and are not intended for  
use with any specific machine or unit.  
This document is provided for customer  
information only, and will not be updated or  
exchanged without customer request.



 **Manufacturer**  
**Dräger Medical Systems, Inc.**  
 3135 Quarry Road  
Telford, PA 18969-1042  
U.S.A.  
 (215) 721-5400  
(800) 4DRAGER  
(800 437-2437)  
**FAX** (215) 723-5935  
 <http://www.draeger.com>

 **EC REP** In Europe, Middle East, Africa, Latin  
America, Asia Pacific distributed by

 **Drägerwerk AG & Co. KGaA**  
Moislinger Allee 53 – 55  
D-23542 Lübeck  
Germany  
+49 451 8 82-0  
**FAX** +49 451 8 82-20 80  
 <http://www.draeger.com>

MS34045 – RI 04 en  
© Drägerwerk AG & Co. KGaA  
Edition: 4 – 2020-05  
(Edition: 1 – 2017-06)



Dräger reserves the right to make modifications  
to the equipment without prior notice