

ICU Medical Plum 360 Documentation

Interfaces: Ethernet

Type: Large Volume Infusion Pump

Manual: See attached Plum360_SOM.pdf

Nmap:

Starting Nmap 7.80 (<https://nmap.org>) at 2023-07-10 16:14 Pacific Daylight Time

Nmap scan report for 192.168.29.177

Host is up (0.0059s latency).

All 65535 scanned ports on 192.168.29.177 are closed

MAC Address: 00:03:B1:45:4E:D8 (Hospira)

Nmap done: 1 IP address (1 host up) scanned in 59.02 seconds

Vulnerability Disclosure Page and/or Contact Email: <https://www.icumed.com/solutions/infusion-therapy/enhance-cybersecurity>; productsecurity@icumed.com